

High-Throughput, Dual-Stream UHPLC/MS/MS Bioanalysis and Data- Deconvolution for Rapid Drug Discovery Applications

62nd ASMS Conference on Mass Spectrometry and Allied
Topics, Baltimore MD 16-June-2014

Brendon Kapinos
Pfizer

1

62nd ASMS

High-Throughput, Dual-Stream UHPLC/MS/MS Bioanalysis and Data-Deconvolution for Rapid Drug Discovery Applications

HDO Group (Hit Discovery and Optimization)

ADME Screening (96- and 384-well assays)

- Metabolism
- DDI & Toxicity
- Disposition
- Permeability & Transporters


In Vitro Bioanalysis

- Transporter Studies
- Enzymology

Technology Development and Validation

2

Our Bioanalytical Work Spans a Broad Scope



3

62nd ASMS

Leveraging Technology to Address Portfolio Needs




Dual-arm (dual-stream) liquid-handler-type LC/MS/MS autosampler

High-throughput, large deck capacity, rapid sample analysis and *flexibility*

4 injection ports per arm, 4 front diverter valves allow multitude of LC modes and/or chemistries, functionalities


4

High-Throughput: "Trap-and-Elute"




62nd ASMS

Continuous Improvement of High-Throughput Bioanalysis




-Optimize Technologies 1.5x5mm trap cartridge
-1.2ml/min LC flow rate (standard trap-and-elute)



-Analytical Sales and Services Sprite Echelon 20x2.1mm
-1ml/min LC flow rate
-Weaker eluent ("Hybrid" trap-and-elute)

"Hybrid" trap-and-elute increases scope of retained compounds and provides separation without impacting throughput

Our Bioanalytical Work Spans a Broad Scope



7


62nd ASMS

Case Study: Transporter Assay Bioanalysis

Challenge	Details
High sample burden	~1200 samples/week (32 compounds x 3 replicates x 6 plates x 2 studies)
Casseted approach	4 compounds/sample (analyte, inhibitor, internal standard, assay control) different polarities, response, LC conditions, <u>requires separation</u>
Diversity of analytes	Ionization efficiency varies
Co-dosed inhibitors	Extremely high concentrations (10-100µM)
Membrane control	Atenolol: small, polar, extreme difference in response across plates
Rapid turnaround	Traditional gradient HPLC: 2min/injection (previously 20 hour run-time)

8


Compound-Dependent Ionization Efficiency



9


62nd ASMS

Compound-Dependent Ionization Efficiency



10


Case Study: Transporter Assay Bioanalysis



11


62nd ASMS

“Dual-stream” autosampling



12


Higher-Throughput Gradient Separation Via Dual-Stream UHPLC



13


62nd ASMS

Higher-Throughput Gradient Separation Via Dual-Stream UHPLC



14


Higher-Throughput Gradient Separation Via Dual-Stream UHPLC



15

62nd ASMS

Higher-Throughput Gradient Separation Via Dual-Stream UHPLC



16

“Dual-Stream” Autosampling

Dual stream analysis is robust and greatly improves throughput

However, column variation over time may lead to signal discrepancy between streams



17


62nd ASMS

Leveraging Flexibility to Expedite Analysis


Solution: Maximize throughput and reproducibility via dual-stream analysis of 2 different analytes simultaneously

Stream 1 injects propranolol

Stream 2 injects diclofenac




Dual-Stream Metabolic Stability (HLM)




62nd ASMS

Dual-Stream Metabolic Stability (HLM)



Dual Stream Analysis: Data Deconvolution




Sound Review (from Sound Analytics)

- ADDA creates batch data file (compound and timing info) during acquisition, file read into Sound Review
- Software only displays analyte collected per stream based on timing data

21


62nd ASMS

Dual Stream Analysis: Data Deconvolution



22

ADDA: Multiple Autosamplers in One



4 injection ports/stream
4 front diverter valves

Flexibility

Multiple LC modes (isocratic, gradient)
Multiple chemistries per LC mode

62nd ASMS

Acknowledgements

Pfizer

Mary Piotrowski
Hui Zhang
Veronica Zelesky
Jefry Shields
Jianhua Liu
Carrie Funk
John Janiszewski
Matt Troutman
Mark West
Bo Feng

Sound Analytics

Wayne Lootsma
Steve Ainley
Will Schramm
Kevin Shirey

Apricot Designs

Jeremy Lawton
Bob Biagioni

AB SCIEX

Tony Romanelli
Joe Doktorski